

BCI practical course : Advanced stimulus presentation

Jason Farquhar (coordinator)
Loukianos Spyrou (Teaching Assistant)

<J.Farquhar@donders.ru.nl>

<l.spyrou@donders.ru.nl>

Learning Goals : Advanced Stimulus Presentation

- Understand:
 - why precise stimulus timing requirements mean we need to run the stimulus generation in a separate process and use more accurate stimulus presentation methods
- Know how to:
 - Use **Psychtoolbox** for more precise auditory/visual stimulus generation and timing

Today's Plan

- Review : Visual Speller Example. Solutions and discussion of problems
- Hands-on : Psychtoolbox based visual speller stimulus

Break

- Test run: 'DCC visitation presentation'

break

- Discussion: Possible projects

Hands on : Visual Matrix Speller

Experiment Task:

- Build a complete visual matrix speller based BCI experiment consisting of 3 blocks:
 - 1) Training/Calibration Block : where the user is presented with matrix speller stimuli and an instruction on which target to attend to
 - 2) Classifier Training Blocks : where the saved labelled data from the calibration block is used to train an ERP classifier
 - 3) Testing Block : where the trained classifier is used predict which symbol the user is attending to and at the end of the sequence this prediction is used to generate feedback

Discussion – Timing

- Timing issues due to:
 - Non-interruptable/single-thread nature of MATLAB
 - **Slow** matlab based drawing command
 - **Variable length** matlab based drawing commands
 - Lack of information about exactly **when** stimulus was presented (just start/end time of `drawnow` command)
 - Once function starts executing it blocks CPU until it's finished..

Solution:

- Use something designed for high speed, high precision stimulus presentation, e.g. **Vision Egg** (Python) or **PyGame** (Python) or **Presentation** etc.
- **Psychtoolbox** – a MATLAB/Octave (mex/OpenGL/PortAudio)

Hands on: Psychtoolbox

- Matlab's basic drawing commands are, both
 - slow (~ 70 ms to execute) and
 - jittery (± 30 ms execution time)
- Psychtoolbox (see www.psychtoolbox.org) allows to direct control of the video/audio hardware
 - Very fast (GPU accelerated) drawing (< 2 ms)
 - Very precise – locked to video hardware (< 2 ms)

Hands on : Psychtoolbox (PTB)

Task:

- Re-write your visual speller BCI from last week to use Psychtoolbox drawing functions
- Compare the timing performance of this version with the Matlab based version

Notes: System architecture

- Again we have a decision:
 1. Keep the drawing code in the experimental control process
 - **Adv**: simple to implement
 - **Dis**: less precise, lots other matlab happening
 2. Move the drawing into a new stimulus generation process
 - **Adv**: more precise timing (less other stuff going on, just draw & send events)
 - **Dis**: more complex

Key Functions : Psychtoolbox

- Note you can get comprehensive help on using PTB functions from there documentation website
 - <http://docs.psychtoolbox.org/Psychtoolbox>
- Or, by using normal matlab help
 - `help Snd`
- Or, by using the functions builtin help, adding a '?' to the command
 - `Screen flip?`

Key Functions : Psychtoolbox

- Screen
 - All PTB drawing commands use this function with it's first argument specifying what exactly to do
 - Screen gives a list of subfunctions,
 - Screen subfunction? Gives help on subfunction
- Open a PTB window at pos with bgColor background
 - `wPtr= Screen('OpenWindow',num,bgColor,rect)`
- Update the display after drawing finished:
 - `screen('Flip',wPtr,when);`
 - If when==0 then wait until next refresh, otherwise wait until when is the system time
 - N.B. Screen **only** changes after a Flip command
- Texture creation and manipulation commands... see next.
- N.B. Add the Psychtoolbox functions to your MATLAB path using:
 - `run ../utilities/initPTBPaths`

Key Concepts : textures

Pyglet uses hardware based textures for very fast drawing.

- A texture is an image which is **pre-loaded** onto the graphics hardware using:
 - `textureID=Screen('MakeTexture',wPtr,image)`
 - `wPtr` is handle to the PTB drawing window
 - `image` is a [w x h x 3/4] image matrix
 - `texelID` is a handle to the created texture
- A texture can then be drawn to the screen **rapidly** using:
 - `Screen('DrawTextures',wPtr,texelID,srcR,destR,angle,filter,alpha,color)`
 - `srcR/destR` specify rectangles [left top right bottom] (ltbr) of the source image to draw at dest position on the screen
- Importantly: at drawing time the texture can be **manipulated extremely rapidly** using the other arguments of 'DrawTextures', e.g.
 - position/size, rotation angle, color, transparency, etc.
- N.B. The screen is only **actually changed** when `Flip` is called!

useful function: mkTextureGrid

PTB equivalent of initGrid to make a grid of textures:

```
[texs,srcRs,destRs]=mkTextureGrid(wPtr,syms)
```

- `wPtr` — PTB window handle (as returned by `Screen('OpenWindow',...)`)
- `syms` - cell array of images or text to layout. The shape of `Syms` gives the shape of the grid
- `texs` - handles to the textures
- `srcRs/destRs` source/destination rectangles (used in `DrawTextures`)
- N.B. Draw the textures using `DrawTextures`, e.g.:

```
Screen('DrawTextures',wPtr,texs,srcRs,destRs,[],[],[],  
[255;255;255]*[0 0 0 1 1 1]);
```

- Makes the first 3 textures **black** (with color [0 0 0]) and the last 3 white (with color [255 255 255]) by changing the color **scaling**, i.e. `Drawn color= originalColor .* drawColorScale`
- Note: in PTB all colors are in `[rgb[a]]` with a range from **0-255**

Summary

- Using PTB (or other dedicated stimulus generation software) means we can get much better stimulus timing accuracy
- BUT Matlab based is simpler, less code / dependencies
 - And for many, e.g. Movement based, improved timing accuracy has no benefit
- So: think about what you **need** before you start!