# BCI practical course: "Hello World" & ERP Viewer

Jason Farquhar (coordinator)

<J.Farquhar@donders.ru.nl>

Loukianos Spyrou (Teaching Assistant)

<l.spyrou@donders.ru.nl>



# Learning Goals

#### **Understand:**

- What is needed to make a BCI, i.e. progress tracking, data acquisition, annotation and processing, stimulus presentation, and an overall process scheduler/sequencer
- How to used event-driven programming ideas coupled to a global shared event pool (blackboard) to provide these facilities
- How the fieldtrip buffer provides the event blackboard which is used for interprocess communication.

#### Know how to:

- What the struct of an 'event' is and how to use it to annotate data with experiment relevant event information
- present simple visual stimulus/feedback to the user/experimenter
- How to wait for specific events, get the necessary data, process it and post the updated results back to the event blackboard
- Startup the buffer and an experiment control Matlab, and how to connect these processes to provide a basic BCI
- Test your experiment with simulated data generated by the signal-proxy
- Debug your experiment when it fails!



# Today's Plan

- Discussion: What do we need to make a BCI?
- Introduction to the Buffer-BCI framework

### break

- Hands-on 1: Hello World
- Hands-on 2: Sequenced Sentences

## break

- Hands-on 3: Visual ERP Viewer
- Brain-test: Visual ERP Viewer




# Discussion: What do we need to make a BCI?

 Based on your prior knowledge and experience with the hands on demo we've did last time.


Discuss: What do we require to make a BCI system?

- Think about:
  - Hardware requirements?
  - Software requirements?
  - Information flows?

## **BCI** information flow




# Gross structure of a typical BCI experiment



# BCI terminology (our group!)

- Epoch/Trial
  - Single BCI prediction
  - e.g. 1-imagined movement, 1 visual-speller flash
- Sequence
  - Short group of epochs (~1min)
  - v. short breaks 1-2sec between epochs (usually automatic)
  - short (usually self-paced) subject break between sequences (~10sec)
- Block/Run/Phase
  - Short group of sequences (>10min)
  - long (~1-2min) subject break between blocks
  - e.g. cap-fitting, calibration, classifier training, on-line use
- Session during one cap-fitting
- Experiment imagined movement, visual-speller etc.

# Flow chart of an individual epoch in a simple BCI experiment



## Requirements: what do you need to build a BCI?

- 1) Way of **tracking** where we are in execution of the experiment flowchart, i.e. block, sequence, epoch number.
- 2)Way of **annotating** data to what the subject was experiencing/doing at that time with what was measured from their brain/body, e.g. LH movement, reading instruction, watching queue, etc.
- **3)**Data acquisition: Drivers to extract data from hardware (and combine data fro different hardware sources)
- **4)**Stimulus Generation: makes stimuli that the subject will experience, for subject instruction, feedback, event-related stimuli
- 5)Something to **process** the signals, firstly to train the classifier, and secondly to decode the users mental state, i.e. do the BCI bit ;-)
- 6) Scheduler (sequencer?) to tie it all these bits together,
  - so the correct functions, i.e. stimulus display, signal processing, are executed
  - at the correct position in the experiment flowchart
  - based on the right bits of measured datacc



# Summary

- To build a BCI we need a system to; track our progress through the experiment, acquire, annotate and process data, present the stimuli and schedule all these processes in an appropriate way.
- Next we introduce a Matlab based system which provides these facilities.

## **Buffer-BCI Framework**

- •We can break the requirements into 4 largely independent communicating processes:
  - 1)Data-acquisitation & annotation
    - ·Get data from hardware
    - Attach annotations (markers, events) to particular data sample
  - 2)Experiment control (scheduling)
    - · Control the flow of the experiment
  - 3)Stimulus generation
    - · Make stimuli when requested by the expt controller
    - · Make feedback based on predictions generated by the sig-processor
  - 4)Signal Processing
    - Process the data based on the annotationss, and generate predictions

## Buffer-BCI Framework

#### Basic Idea:


- set of independent processes
- any process can send/recieve dataannotation events

events are visible to all other processes

 Processes communication implemented by sending recieving events

 (N.B. As all events are saved with the data, annotations are automatically archived for later off-line use).

 Similar in concept to that used in 'Blackboard architectures' for AI, see < en.wikipedia.org/wiki/Blackboard\_syste



Computation

# Ft-buffer based Implementation

- Buffer-BCI framework implemented using the fieldtrip-buffer system (fieldtrip.fcdonders.nl/development/realtime)
- Ft-buffer provides:
  - Drivers for data-acquisation
  - 1)buffer storage for data (~last 1 minute data)
  - 2) buffer storage for events (~last 50 events)

#### Idea:

- Buffer events store represents the blackboard used for inter-process communication (IPC)
- Every event has timestamp (sample number) used for data-annotation

# (Key concept) event structure

#### sample

 time at which event occurred in samples from start of experiment

#### type

arbitrary event type (usually a string)

#### value

 arbitrary event value (usually string or number)

#### duration (optional)

duration of the event in samples

#### offset (optional)

- zero-time for the event.
- Usually, offset from sample at which the event actually started.

```
Examples:
Visual speller "flash";
ev=struct('sample',123,...
'type','stimulus.flash',...
'value',[0 0 1 0 0],...
```

'offset',0,'duration',0)

```
Classifier prediction:
ev=struct('sample',123,...
'type','prediction',...
'value',[-1 -1 -1 1 -1],...
'offset'.0.'duration',0)
```

```
Imagined Movement event:
ev=struct('sample',123,...
'type','stimulus.move',...
'value','left-hand',...
'offset',0,'duration',300)
```

Compact notation: s:123,t:'stimulus.flash',v:[0 0 1 0 0],o:0,d:0

# (key functions) Event manipulation

evt=mkEvent(type,value,[sample,offset,duration])

make a buffer event, with sensible defaults

#### sendEvent

- evt=sendEvent(type,value,[sample,offset,duration,host,port])
- evt=sendEvent(evt,[host,port])
- Send event to the buffer on machine host at port.

#### *mi*=matchEvents(*evts*, *mtype*, *mval*)

- Find events with type mtype and value mval in evts a vector of event structures.
- mtype can be cell-array of types to match, e.g. {'type1' 'type2'}
- mval can be cell-array of values to match, e.g. {'val1' 10 'val3'}
- match if any mtype matches and any mval matches,
  - i.e. above matches (t:'type1',v:10), (t:'type2',v:10),(t:'type1',v:'val1')
- mi is logical vector of which evts matched
- N.B. Empty ([]) or '\*' mtype/mvalue matches everything





## Hands-on 1: "Hello World"

## **Experiment Task**

- Display the string "Hello World" (or any other prespecified string) on the screen, and wait for a key to be pressed to exit
- Send events to annotate what has happened, e.g. startup, string display, key-pressed, shutdown etc.

### Method:

- Start from the 'helloworld-skel.m' function skeleton
  - contains initialisation code to connect to the ft\_buffer
  - Some examples of functions you may find useful



# Running buffer-bci code:

- You need to have (at least) the following processes running:
  - 1) fieldtrip-buffer and data-acquisation:
    - simulated data: separate buffer and simulated acquisation:
      - dataAcq/startBuffer.sh
      - dataAcq/startSignalProxy.sh
    - Biosemi : combined buffer and data-acquisation:
      - dataAcq/startBiosemi.sh
    - Emotiv: combined buffer and data-acquisation:
      - dataAcq/startEmotiv.sh
    - Mobita: separate buffer and data-acquisation
      - dataAcq/startBuffer.sh
      - Virtual Windows machine for the data-acquisttion client (Polybench)
  - 2) experiment control process



# Running buffer-bci code:

- You need to have (at least) the following processes running:
  - 1) fieldtrip-buffer and data-acquisation:
  - 2) experiment control process
    - Usually your .m file in matlab
    - (N.B. May also contain the stimulus-presentation and signal-processing code – as in this exercise)
    - Start MATLAB
    - connect to buffer (see header in helloworld-skel.m)
    - run control process

# Useful (debugging) Functions:

 Seeing the what events are sent and when is important for debugging experiment flow.

[]=eventViewer(host,port,mtype,mval)

- Print all events matching (mtype, mval)
- (mtype,mval) are as used in matchEvents

```
N.B. you can run eventViewer directly using:
utilities/startEventViewer.{bat,sh}
```

## Note: event timestamps

- Accurate event time-stamps are critical for evoked potential analysis
  - >10ms event jitter causes significant reduction in signal quality
- However,
  - data-acquisation may only send data every >20ms
  - And this data may be subject to additional network delays of >20ms
- Stop this jitter reducing time-stamp accuracy by;
  - aligning (and tracking) computers real-time-clock and datasample clock to prevent this jitter reducing time-stamp accuracy

## **Hands-on 2: Sequenced Sentences**

#### **Experiment Task**

- display set of sentences on the screen where every second 1 more character gets added to the sentence
- pause for 5 seconds between sentences (and/or wait for key press)
- send events for everything that happens

#### Assignment:


- Make flowchart
- Write code -> test -> debug -> until it works :-)
- Start from runSentences-skel.m

#### **Useful Functions:**


- sleepSec(time)
  - cause matlab to sleep for the indicated duration in seconds (more accurate than 'pause').



# Flowchart: sequenced sentences



# Events and processing functions





## Hands-on 3: Event Echo-Server

## **Events for IPC**

- As well as being used for data annotation, Events are used for inter-processcommunication,
  - e.g. to communicate results from signal-processing to stimulus presentation
- To use events in this way, each process must
  - monitor for new events
  - filter out the events it should react to
  - send response events



## **Hands-on 3: Event Echo-Server**

## **Experiment Task**

- Write a simple echo-server which:
  - Waits for any incomming event, and
  - Responds by sending a 'echo' event with the same value but type='echo'
  - Quits if it recieves and event with type='exit'
- N.B. Don't 'echo' your own echo events!

## **Assignment:**

start from : echoServer-skel.m

N.B. send 'keyboard' events by pressing keys in the signal proxy window.



## **Useful Functions:**

[devents,state]=...
buffer\_newevents(host,port,state,mtype,mval,timeout\_ms)

- wait for any new events matching (mtype, mval)
  - Matching done by matchEvents

```
mtype – can be cell-array of types to match, e.g. {'type1' 'type2'} mval – can be cell-array of values to match, e.g. {'val1' 10 'val3'} match if any mtype matches and any mval matches
```

- return the matched events in the vector of structure(s) devents
- state is the match state, used to track which events have been processed between function calls
- Return after timeout\_ms milliseconds even if no matching events found

# Echo-Server in different languages

- Basic echo-server example has been implemented in multiple languages
- example directory contains example implementations in different languages
- MATLAB: matlabclient/matlabclient.m
- JAVA: javaclient/javaclient.java
- C#:csharpclient/csharpclient.java
- Python: pythonclient/pythonclient.py
- C:cclient/cclient.c



## Hands-on 3: ERP Viewer

#### **Experiment Task**

- In 5 sequences of 10 seconds:
  - Every 1 seconds: either randomly display or don't display a cross (+) on the screen for 200ms
- Display a 'Press key to continue string' between sequences, and wait for key press to move to the next sequence
- For every 'stimulus event', i.e. point when the '+' could have been displayed, record 600ms of data annotated with whether it was a '+' or not
- Every time you get some data, compute an average of the EEG data for that type of stimulus, i.e. + or no-+, and display the resulting averages as a multi-plot on the screen
- N.B. You will need a separate signal processing process to: get the data, compute the ERP and display the results!

#### Assignment:

- Make flowchart
- For the expt-control & simulus presentation start from : runStimulus-skel.m
- For the signalProcessing & results genration use: runSigProc-skel.m



## **Useful Functions:**

[data,devents,state]=buffer\_waitData(host,port,state,...
'startSet',startEvts,'trlen\_samp',samp,'exitSet',exitEvts)

- for all events matching startEvts record samp samples of data
- until an event matching exitEvts is generated
- *startEvts* and *exitEvts* specify the events to match in the format:
  - type event type has this value
  - {'type' val} event has type==type and value==val
  - {{'type1' "type2'}} event has type == 'type1' or 'type2'
  - {{'type1' 'type2'} {val1 val2}} event has type == 'type1' or 'type2' and value== val1 or val12
- return the matched event structure(s) in *devents* and corrospending data in *data* 
  - Data is a vector of structures. data.buf = [nChannels x nSamples] raw EEG data
- state is the match state, used to identify which events have been processed between function calls
- N.B. ExitEvts has the special event type 'data' which returns as soon as the data is available for the first matched startEvt



# Summary

- BCI can be broken into 4 processes: dataacquisation, experimental control, signal processing, and stimulus presentation
- buffer\_bci framework : uses buffer events as a blackboard for inter-process communication